

TOPIMPLANT



Product Catalog

产品手册 -3D 替代体

Aidite DentalGroup

The background features a light gray grid of circles that curves and fades into the background. In the lower half, there are several overlapping, wavy, translucent bands in shades of gray and white, creating a sense of depth and movement.

让种植修复更简单

Contents

产品目录

公司介绍 01

关于我们 01

公司使命愿景 03

产品介绍 05

3D替代体 05

桥体3D替代体 17

数字替代体安装夹具 21

About Us

关于我们

极植科技(深圳)有限公司成立于2022年,注册资本1500万元,公司位于广东省深圳市宝安区,是一家集科研、开发制造、销售于一体的牙科种植配件及CAD/CAM数字化解决方案提供商。

公司前身是深圳市三专医疗科技有限公司数控加工中心,专注于种植配件加工领域8年的时间,现有威力铭、西铁城等20多台加工及检测设备,主要产品有替代体、扫描杆、转移杆、技工螺丝、定位环、螺丝起、扭力扳手及各类种植辅助工具。公司视产品质量为生命,重点打造了自动化检测流水线,保证产品外观及螺丝100%检验,确保了产品高质量和稳定性。公司拥有在数控加工领域20多年经验的资深专家,坚持产学研合作。

2015

成立数控加工中心,购置瑞士威力铭五轴五联动加工中心,成为第一批使用高精尖威力铭机床切削种植配件的企业。

2017

购置第一批西铁城走芯机,成为较早同时拥有瑞士威力铭和日本西铁城设备的种植配件加工企业。



2020

2020年公司完成了覆盖欧系、韩系、美系等60多个种植系统1000多种型号种植配件的开发工作。

2022

2022年从三专公司独立出来成立极植科技(深圳)有限公司,全面启动基台柱、个性化基台、个性化复合基台、个性化基台桥、个性化基台杆卡等三类医疗器械的注册工作。

2024

我们取得了个性化基台及配件,可切削基台柱在中国的注册证。

Company Mission And Vision

公司使命愿景

愿景 Vision

成为世界顶级的牙科种植配件加工中心

使命 Mission

让种植修复更简单

核心价值观 Core Values

为客户创造价值

以奋斗者为本

自省 创新 诚信 合作





Digital Analog

3D替代体

精心设计

- 其独特的双抗旋设计、更小的外形以及更短的长度,都体现了我们对细节的考虑和精心设计

就位保障


- 底部的螺丝结构确保了替代体纵向就位,增强了稳定性和精度

反馈设计

- 其独特的凹槽设计使得当替代体插入打印的数字模型时,用户可以通过清晰的“咔哒”声音反馈感知其是否准确安装

易识别

- 替代体上的颜色代码和激光雕刻标识帮助医生快速、准确地识别不同的替代体,避免技工选择错误的替代体
- 这些特点展现了我们数字替代体的高质量和用户友好的设计,旨在提供一个准确、方便和安全的替代体



Digital Analog

3D替代体

种植品牌	规格型号	平台大小	产品编码	数据库编码
Dentium	SuperLine Implantium II	■ 3.4	AD DEN	TP01
		■ Regular	AD DER	
	NR Line	■ 3.1/3.6/4.3/5.0	AD DNR	
	Simple Line	■ 4.8	AD DSR	
		■ 6.5	AD DSW	
OSSTEM / Hiossen	TS/ET	■ Mini	AD OTM	TP02
		■ Regular	AD OTR	
	SS	■ Regular	AD OSR	
		■ Wide	AD OSW	
	US	■ Mini	AD OUM	
		■ Regular	AD OUR	
		■ Wide	AD OUW	
MegaGEN	EZ Plus	■ Mini	AD MEM	TP03
		■ Regular/Wide	AD MER	
	MiNi	■ 3.0	AD MM3	
	AnyOne	■ Narrow	AD AON	
		■ Regular	AD AOR	
	AnyRidge	■ 3.5/4.0/4.25/4.5/5.0	AD MAR	
	ST	■ Mini	AD MSS	
		■ Standard	AD MSM	
Cortex	Internal Hex	■ Regular	AD CTX	TP04
	Conical Platform	■ NP	AD CCN	
		■ RP	AD CCR	
		■ WP	AD CCW	

Digital Analog

3D替代体

种植品牌	规格型号	平台大小	产品编码	数据库编码
BEGO	S/RI/PS/SB/SBX/RS/R SX/RI	■ 3.25/3.75	AD BGM	TP05
		■ 4.1	AD BGS	
		■ 4.5	AD BGL	
		■ 5.5	AD BGW	
DIO	SM(TORX)	■ Narrow	AD DSN	TP06
		■ Regular/Wide	AD DSR	
	UF(II)	■ Narrow	AD DUN	
		■ Regular/Wide	AD DUR	
Straumann	Bone Level	■ SC	AD SBS	TP07
		■ NC	AD SBN	
		■ RC	AD SBR	
	Tissue Level	■ NNC	AD SSC	
		■ NN	AD SSN	
		■ RN	AD SSR	
		■ WN	AD SSW	
	BLX	■ RB/WB	AD BLX	
Bicon	/	■ 2.0	AD BC20	TP08
		■ 2.5	AD BC25	
		■ 3.0	AD BC30	
DENTIS	S-Clean	■ Narrow	AD DTS N	TP09
		■ Mini	AD DTS M	
		■ Regular	AD DTS R	

种植品牌	规格型号	平台大小	产品编码	数据库编码
Zimmer	TSV	3.5	AD ZTG	TP10
		4.5	AD ZTP	
		5.7	AD ZTY	
	Tapered Swiss Plus	4.8	AD ZRN	
	Spline	3.25	AD ZSN	
		3.75/4.0	AD ZSR	
		5.0	AD ZSW	
	Azure C-System	S	AD AZS	
		L	AD AZL	
MIS	SEVEN	NP	AD M7N	TP11
		SP	AD M7S	
		WP	AD M7W	
	C1	NP	AD C1N	
		SP	AD C1S	
		WP	AD C1W	
Anthogyr	Axiom BL	BL	AD AOM	TP12
		2.8	AD AOM28	
Nobel Biocare	Replace	NP	AD NRN	TP13
		RP	AD NRR	
		WP	AD NRW	
	Active	3.0	AD NA3	
		NP	AD NAN	

Digital Analog

3D替代体

种植品牌	规格型号	平台大小	产品编码	数据库编码
Nobel Biocare	Active	 RP	AD NAR	TP13
		 WP	AD NAW	
	Branemark	 NP	AD NBN	
		 RP	AD NBR	
		 WP	AD NBW	
Neoss	ProActive	 3.25	AD NSN	TP14
		 3.5	AD NSR	
Dentsply	Astra OsseoSpeed TX	 3.0	AD DA3	TP15
		 Narrow(3.5/4.0)	AD DAS	
		 Regular(4.5/5.0)	AD DAL	
	Astra EV	 3.0	AD AEG	
		 3.6	AD AEP	
		 4.2	AD AEY	
		 4.8	AD AEB	
		 5.4	AD AE5	
	Ankylos	 C/	AD ANK C/X	TP16
		 /X		
	Xive	 3.0	AD DX3	TP17
		 3.4	AD DXG	
		 3.8	AD DXY	
		 4.5	AD DXB	
		 5.5	AD DXR	

种植品牌	规格型号	平台大小	产品编码	数据库编码
CAMLOG	Camlog	■ 3.3	AD CCG	TP18
		■ 3.8	AD CCY	
		■ 4.3	AD CCP	
		■ 5.0	AD CCB	
		■ 6.0	AD CC6	
	isy	■ 3.8/4.4/5.0	AD ISY	
CAMLOG	Conelog	■ 3.3	AD CON33	TP19
		■ 3.8/4.3	AD CON38/43	
		■ 5.0	AD CON50	
BIOMET 3i	Certain	■ 3.4	AD BCP	TP20
		■ 4.1	AD BCB	
		■ 5.0	AD BCY	
		■ 6.0	AD BCG	
Medentis Medical	ICX	■ 3.3	AD ICX33	TP21
		■ Regular	AD MIC	
SIC	Internal Hex	■ 3.3	AD SHB	TP22
		■ 4.2	AD SHP	
Neobiotech	ISIS/IS II/IS-II Active	■ Narrow	AD NEN	TP23
		■ Regular/Wide	AD NER	
		■ Wide	AD NEW	
B&B DENTAL	Internal Hex	■ Mini	AD BBM	TP24
		■ Regular	AD BBR	

Digital Analog

3D替代体

种植品牌	规格型号	平台大小	产品编码	数据库编码
C-TECH	Narrow/EstheticLine	■ EL	AD CTE	TP25
		■ ND	AD CTN	
Alpha-BioTec	IH	■ 3.3/3.75/4.2/5.0	AD ABT	TP26
	CHC	■ 3.2/3.5	AD CHC	
	CS	■ 3.75/4.2/5.0	AD CS	
Biodenta	Bone Level	■ B0	AD BDT-B0	TP27
		■ B1	AD BDT-B1	
		■ B2	AD BDT-B2	
ImplantDirect	InterActive	■ 3.0	AD IDA30	TP28
		■ 3.4	AD IDA34	
	Legacy	■ 3.0	AD IDL30	
		■ 3.5	AD IDL35	
		■ 4.5	AD IDL45	
		■ 5.7	AD IDL57	
AB-dent	Two Piece	■ Standard	AD ABS	TP29
		■ Narrow	AD ABN	
Adin	RS	■ Regular	AD ARS	TP30
	CloseFit	■ UNP	AD ACU	
		■ NP	AD ACN	
		■ RP	ADACR	
		■ WP	AD ACW	

种植品牌	规格型号	平台大小	产品编码	数据库编码
BioHorizons	Internal Hex	 3.0	AD BIOH-IH30	TP31
		 3.5	AD BIOH-IH35	
		 4.5	AD BIOH-IH45	
		 5.7	AD BIOH-IH57	
	External Hex	 3.7	AD BIOH-EH37	
		 4.2	AD BIOH-EH42	
		 5.2	AD BIOH-EH52	
		 6.2	AD BIOH-EH62	
CLC	CONIC	 Regular	AD CCR	TP32
Bredent	SKY	 3.5N/4.0/4.5/5.5	AD SKY	TP33
Ritter	Snap SB/LA	 Narrow	AD RSN	TP34
		 Standard	AD RSS	
SG	SST/SRT	 Narrow	AD SGN	TP35
		 Regular	AD SGR	
		 Wide	AD SGW	
		 Wide6.0	AD SG6	
Ticare(MG)	Inhex Quattro	 MINI	AD MGM	TP36
		 STD	AD MGR	
		 MAXI	AD MGW	
Keystone	Genesis	 3.5/3.8	AD KGM	TP37
		 4.5	AD KGR	

Digital Analog

3D替代体

种植品牌	规格型号	平台大小	产品编码	数据库编码
Keystone	Genesis	■ 5.5/6.5	AD KGW	TP37
	Prima(TiLobe)	■ SD	AD KPS	
		■ RD	AD KPR	
		■ WD	AD KPW	
Thommen	SPI	■ 3.5	AD TSY	TP38
		■ 4.0	AD TSG	
		■ 4.5	AD TSB	
		■ 5.0	AD TSW	
		■ 6.0	AD TSP	
Mytis Arrow	One-Stage	■ Normal	AD MAN	TP39
		■ Wide	AD MAW	
HumanTech	RatioPlant	■ Mini	AD HTM	TP41
		■ Regular	AD HTR	
		■ Wide	AD HTW	
TBR	Z1	■ Mini	AD TZM	TP42
		■ Infinity/8	AD TZ8	
	Bone Level	■ Mini	AD TBM	
		■ Infinity/8	AD TB8	
Anker	Bone Level	■ Narrow	AD ABN	TP43
		■ Standard	AD ABS	
	Tissue Level	■ NN	AD ATN	
		■ RN	AD ATR	

种植品牌	规格型号	平台大小	产品编码	数据库编码
Anker	Tissue Level	■ WN	AD ATW	TP43
Neodent	GM	■ 3.5/3.75/4.0/4.3/5.0/6.0	AD NGM	TP44
	CM	■ 3.5	AD NCM	
		■ 4.3		
		■ 5.0		
WEGO	/	■ Mini	AD WGM	TP45
		■ Regular	AD WGR	
JD	JDEvolution Plus	■ JDE	AD JDE	TP46
BlueSkyBio	MAX	■ NP	AD BMN	TP47
Kyocera	Bone Level	■ NP	AD KBN	TP49
		■ RP	AD KBR	
		■ WP	AD KBW	
Global D	In-Kone	■ Regular	AD GDI	TP50
EBI	BLT	■ Mini	AD EBM	TP51
		■ Regular	AD EBR	
	Tissue Level/Precision/Octa	■ Internal Octa4.8	AD EBT	
DYNA	Helix	■ 3.2	AD DHW	TP52
		■ 3.6	AD DHG	
		■ 4.2	AD DHB	
		■ 5.0	AD DHY	
LASAK	BioniQ/BioniQ Plus	■ 2.9	AD LBY	TP53
		■ 3.5	AD LBB	

Digital Analog

3D替代体

种植品牌	规格型号	平台大小	产品编码	数据库编码
LASAK	BioniQ/BioniQ Plus	■ 4.0	AD LBG	TP53
		■ 5.0	AD LBP	
Paltop	Internal Hex	■ Narrow	AD PIN	TP54
		■ Standard	AD PIS	
		■ Wide	AD PIW	
Shinhung	Luna	■ NP	AD SLN	TP55
		■ RP	AD SLR	
Warantec	IT	■ Regular	AD WIT	TP56
	IU	■ Mini	AD WUN	
		■ Regular	AD WUR	
COWELL INNO	Submerged	■ Mini	AD INM	TP57
		■ Regular	AD INR	
SNUC	AF	■ Mini	AD SNM	TP58
		■ Regular	AD SNR	
ZUGA	Bone Level	■ 3.5/4.3/5.0	AD ZGB	TP59
Biotech Dental	Kontakt	■ 3.0	AD BT3	TP60
		■ 3.6/4.2/4.8/5.4	AD BTH	
	WeegoFR	■ 3.6/4.2/4.8	AD BTWG	
DENTIN	Class/Rapid/Prestige	■ 3.3/3.8/4.2/5.0/6.0	AD DTY	TP61
Biotem	AR	■ Narrow	AD BTN	TP62
		■ Regular	AD BTR	
SGS	P5D/P5G/P1D/P7D	■ Regular	AD SGS R	TP63

种植品牌	规格型号	平台大小	产品编码	数据库编码
Sweden & Martina	Premium	 3.3	AD SM300	TP64
		 3.8	AD SM380	
		 4.25	AD SM425	
		 5.0	AD SM500	
Intralock	Internal Hexagon	 3.5	AD ILK35	TP65
Bio3	ADVANCED	 Standard	AD B3S	TP66
		 Wide	AD B3W	
GMI	Avantgard	 NP	AD GAN	TP67
		 RP	AD GAR	
		 WP	AD GAW	
IML	UNIVERSE	 NP	AD MLN	TP68
		 RP	AD MLR	

Multi-Unit Digital Analog

桥体3D替代体



Multi-Unit Digital Analog

桥体3D替代体

种植品牌	规格型号	产品描述	订货编号	数据库编码
Dentium	SuperLine	Screw Abutment 4.5	MAD DT45	TP01
		Screw Abutment 5.5	MAD DT55	
Osstem	GS/TS	Multi unit	MAD OST	TP02
		Convertible Abutment Mini/Regular Ø4.0	MAD CA40	
		Convertible Abutment Regular Ø5.0	MAD CA50	
		Convertible Abutment Regular Ø6.0	MAD CA60	
MegaGEN	AnyOne/AnyRidge	N-Type	MAD MGG-N	TP03
		S-Type	MAD MGG-S	
Cortex	Internal Hex	Multi-Unit Hex Abutments	MAD CTX	TP04
BEGO	/	Multi Plus	MAD BG	TP05
DIO	UF(II)	4.8	MAD UF48	TP06
Straumann	/	Screw-retained Abutment NC Ø3.5	MAD SRA3.5*	TP07
		Screw-retained Abutment NC Ø3.5 Single Crown	MAD SRA3.5	
		Screw-retained Abutment NC/RC Ø4.6	MAD SRA4.6*	
		Screw-retained Abutment NC/RC Ø4.6 Single Crown	MAD SRA4.6	
		SynOcta RN Abutment	MAD RN	
		SynOcta WN Abutment	MAD WN	
		Multibase NC Ø3.5	MAD MB35	
		Multibase NC/RC Ø4.5	MAD MB45	
		Multibase RC Ø6.5	MAD MB65	
Bicon	/	Fixed-Detachable Abutments	MAD BC	TP08
Dentis	S-Clean	MU Abutment System	MAD DTS	TP09

Multi-Unit Digital Analog

桥体3D替代体

种植品牌	规格型号	产品描述	订货编号	数据库编码
MIS	SEVEN	Ø4.8	MAD MIS	TP11
Anthogyr	Axiom BL/TL	Ø4.0mm MULTI-UNIT RANGE	MAD ANT40	TP12
		Ø4.8mm MULTI-UNIT RANGE	MAD ANT48	
Nobel Biocare	Active/Replace /Branemark	MU-NP/RP	MAD NB	TP13
		MU-NP/RP Plus	MAD NB V2	
		MU-WP	MAD NBW	
Neoss	ProActive	Access Abutment	MAD NES	TP14
Dentsply	Astra TX	Uni Abutment 20°	MAD AT20	TP15
		Uni Abutment 45°	MAD AT45	
	Astra EV	Uni Abutment EV 33°	MAD AE33	
		Multibase	MAD AEB	
	Ankylos	4.2	MAD ANK	TP16
	Xive	TG	MAD DXV	TP17
Camlog	COMFOUR	Ø4.3	MAD COM43	TP18
		Ø6.0	MAD COM60	
BIOMET 3i	Certain	Single and multi-unit screw-retained Restorations	MAD 3i	TP20
SIC	SICace/SICmax/SICtapered	Safe on Four	MAD SIC	TP22
Neobiotech	IS	Ø4.8	MAD NEO	TP23
BB-dent	/	MULTI-USE ABUTMENTS	MAD BB	TP24
Alpha-BioTec	IH	Multi Unit Abutment System	MAD ABT	TP26
AB-dent	/	single unit	MAD AB	TP29
Adin	RS	TMA	MAD ARS	TP30

种植品牌	规格型号	产品描述	订货编号	数据库编码
Ticare(MG)	Inhex	4.0, 10°	MAD MG10	TP36
		4.0, 30°	MAD MG30	
		4.8, 20°	MAD MG20	
Neodent	GM	Ø4.8	MAD NGM	TP44
Kyocera	Bone Level	4.8	MAD KYC	TP49
Warantec	IU	Ø4.0	MAD WR40	TP56
		Ø5.0	MAD WR50	
COWELL	INNO	Ø4.5	MAD CWN45	TP57
		Ø5.5	MADCWN55	

Digital Analog Positioning Jig

数字替代体安装夹具

- 精密制造, 手感舒适
- 兼容多个种植系统
- 颜色分类标记, 一目了然
- 可拧替代体配套的底部锁紧螺丝
- 兼容登腾复合基台螺丝刀
- 高强度不锈钢材料制造, 坚固耐用



M2.0 黄色 				
Astra TX 4.5/5.0	DIO SM R/W	DIO UFII R	MegaGen EZ R	Biodenta B0
Nobel Active RP/WP	Nobel Replace RP/WP	Camlog 5.0/6.0	Neobiotech R	Dentium
SG R	INNO	Straumann RN/WN	Osstem Mini	Osstem TS R/Hiossen ET R
M1.8 蓝色 				
DIO SM Mini	MegaGen Anyridge	Straumann NN	AlphaBioTec IH	MIS
BEGO	Cortex	Adin	BioHorizons	BB R
Zimmer TSV	Nobel Replace NP			
M1.6 紫色 				
Astra TX 3.5/4.0	MegaGen EZ Mini	Nobel Active NP	Straumann NC/RC	Xive 3.4/3.8/4.5
Camlog 3.3/3.8/4.3	3i CER 3.4/4.1/5.0	MIS	SPI 4.0/4.5/5.0	ICX
M1.4 黑色 				
MegaGen Mini	DIO UFII Narrow	Dentis Narrow	Straumann SC	Astra TX 3.0
Astra EV 3.0	Nobel Active 3.0			
M1.4 灰色 (MUA)				
Nobel	Osstem/Hiossen	MegaGen	Straumann SRA4.6	Xive
Cortex	ICX	BioHorizons	Adin	AB
MIS	Zimmer			

TOPIMPLANT

Aidite DentalGroup

☎ (+86) 188-3358-0801 (微信同号)

✉ sales@topimplant.cn | info@topimplant.cn

🌐 www.topimplant.cn

📍 深圳市宝安区福海街道和平社区桥和路306号福海国际科技园A区A1栋138 极植科技



Focus more...